



## Equipment Decontamination

### 1. Purpose:

---

To provide instruction on how to properly and safely decontaminate equipment

### 2. Scope:

---

Applies to everybody working in CL2 Lab

### 3. Prerequisites:

---

WHMIS and Laboratory Biosafety Training (EHS601)

### 4. Responsibilities:

---

It is the responsibility of the PI and lab personnel to ensure that equipment is properly decontaminated, and in the case of transport, properly prepared.

### 5. Personal Protection Equipment (PPE):

---

Lab coat, nitrile gloves, safety goggles/glasses



### 6. Procedure:

---

- 1) Before beginning decontamination:
  - Equipment should be unplugged
  - All materials should be removed and properly discarded or stored
  - Refrigerators should be allowed to reach room temperature prior to decontamination and freezers must be defrosted
  - All materials used to soak up thawed water should be collected and treated as biohazardous waste
  - Incubators and other equipment with water temperature jackets must be drained if they are to be moved and the water drained from the jacket should be treated as biohazardous waste
- 2) In most cases, the following decontamination and contact times are appropriate:
  - *Tissue specimens, viruses (non-retroviral viruses) and non-spore forming bacteria:* 1% sodium hypochlorite for 30 minutes, or 70% ethanol for 20 min (read manual to choose compatible disinfectant)

- *Spore forming bacteria*: 1% sodium hypochlorite (freshly prepared) for  $\geq 30$  min
  - *Fungi*: 1% sodium hypochlorite (freshly prepared) for 20 min
  - *Others*: Virox, Clidox, gaseous Formaldehyde or H<sub>2</sub>O<sub>2</sub> - discuss with the Biosafety Officer
- 3) Prepare a sufficient amount of decontaminant and transfer to a squeeze bottle (do not use spray bottle)



- 4) Squeeze covering all inner surfaces of the item with the decontaminant including door, handles, etc., and leave for the appropriate contact time
- 5) Wipe with absorbent material, collect and treat as biohazardous waste
- 6) Repeat steps 4 and 5 if the surfaces are particularly covered with dirt or “grime”
- 7) Following cleaning with 1% sodium hypochlorite, rinse well with water to prevent corrosion
- 8) Once the equipment is decontaminated, label with a “Decontaminated” sign that includes the date and the name of the person who was responsible for the decontamination process.